

## Hearing Aid Comparison



This hearing aid is best for	1:1 conversations, easy to maintain, self programmable	1:1 conversations, conversations in small groups in noise, hearing the TV	General conversations in background noise, quiet meeting rooms, busy meeting rooms, hearing the TV and phone	General conversations in background noise, quiet meeting rooms, busy meeting rooms, hearing the TV and phone and also rechargeable	Hearing in complex listening environments and background noise, listening to music, TV, phone
Pricing	FROM \$1,900pp	FROM \$1,995pp	FROM \$5,350pp	FROM \$5,650pp	FROM \$6,590pp
Bluetooth	✗	✓		✓	✓
Battery Type	Disposable (312)	Disposable (13)	Disposable (312)	Rechargeable	Rechargeable
Battery Life	7-10 days	7-10 days	7-10 days	18 hours	18 hours
Voice prompts	✓	✗		✗	✗
Dimensions	23 x 10 x 6mm	29 x 12 x 8mm	27 x 11 x 7mm	28 x 12 x 8mm	28 x 12 x 8mm
Water repellent nano-coating	✓	✓		✓	✓
Telecoil	✓	✓		✓	✗
Volume switch	✓	✓		✓	✓
Support App	ihearyou	myphonak		myphonak	myphonak
App Compatibility - iOS	iOS 10.2	iOS 12		iOS 12	iOS 12
App Compatibility - Android	4.1	Android OS 7.0		Android OS 7.0	Android OS 7.0
Direct calls/ audio streaming	✗	Bluetooth® 4.2 wireless technology and most older Bluetooth phones			
Tinnitus Balance	✗	✓		✓	✓