





How to choose the right tube length for Marvel

Grab a ruler and get ready to measure your best speaker tube length. Make sure you measure from 0 where your ear joins to your head.



You can use a pair of glasses
Measure the distance from the underside of your glasses frame (A) to the top third of your ear canal (B)



You can use a sheet of paper
Rest a piece of paper horizontally on top of your ear, and measure from the bottom of the paper (A) to the top third of your ear canal (B)

Marvel hearing aid tube length guide

MEASUREMENT	SIZE
19-23 mm	0
23-27 mm	1
27 - 34 mm	2



If you fit within two ranges, pick the larger size.