

Product Comparison Chart

	Phonak Audéo™ Marvel 30-13T	Phonak Audéo Paradise 70-312	Phonak Audéo Paradise 70-R	Phonak Audéo Paradise 90-R
This hearing aid is best for	1:1 conversations, conversations in small groups in noise, hearing the TV	General conversations in background noise in stationary situations, quiet and busy meeting rooms, hearing the TV and streaming from other Bluetooth enabled smart devices as well as enjoyment of music.	General conversations in background noise in active and stationary situations, quiet and busy meeting rooms, hearing the TV and streaming from other Bluetooth enabled smart devices as well as enjoyment of music	Hearing in the broadest range of active and stationary listening environments including different levels of reverberation, background noise, speech levels, music enjoyment and car situations.
Pricing				
Price per pair**	\$ 1,995.00	\$ 5,880.00	\$ 6,200.00	\$ 6,990.00
Bluetooth				
Bluetooth	Yes	Yes	Yes	Yes
Bluetooth device memory	2	8	8	8
Features & Specifications				
Battery type	Disposable Zinc Air (13)	Disposable Zinc Air (312)	Lithium Ion Rechargeable	Lithium Ion Rechargeable
Battery life	3-4 days	4 days^	24 hours*	24 hours*
Dimensions	29 x 12 x 8mm	27 x 11 x 7mm	29 x 12 x 8mm	29 x 12 x 8mm
Acoustic tube or Receiver	Receiver	Receiver	Receiver	Receiver
Advanced environmental noise management	Yes	Yes	Yes	Yes
Advanced wind noise management	No	Yes	Yes	Yes
Advanced feedback management (to prevent whistling)	Yes	Yes	Yes	Yes
Impulse noise suppression	No	Yes	Yes	Yes
Voice prompts	No	No	No	No
Tap Control / voice assistant	No	No	Yes	Yes
Motion Sensor Hearing	No	No	Yes	Yes
Telecoil	Yes	No	No	No
Tinnitus	Yes	Yes	Yes	Yes
Volume switch	Yes	Yes	Yes	Yes
Direct calls / audio streaming	Bluetooth® 4.2 wireless technology and most older Bluetooth phones			
Water repellent nano-coating	Yes	Yes	Yes	Yes
App compatibility				
Support app	myPhonak	myPhonak	myPhonak	myPhonak
Real time adjustments via app	Yes	Yes	Yes	Yes
App compatibility - iOS	iOS 12	iOS 12	iOS 12	iOS 12
App compatibility - Android	Android OS 7.0	Android OS 7.0	Android OS 7.0	Android OS 7.0
Accessories				
TV streamer compatibility	Yes	Yes	Yes	Yes
Compatible accessories	TV Connector	TV Connector	TV Connector	TV Connector
	PartnerMic	PartnerMic	PartnerMic	PartnerMic
	RemoteControl	RemoteControl	RemoteControl	RemoteControl
**Single devices also available		^Mixed usage including streaming	*Streaming will effect battery life	*Streaming will effect battery life